



DORNIER
CASCADE™
URETERAL STENT

Our hydrophilic coated CASCADE™ stent with double pig tail design and one step insertion helps to facilitate comfort and stability

DRAIN


Activated by fluid, the hydrophilic coating creates a smooth surface between the stent and tissue for easy insertion and passage through delicate ureteral tissue

Designed with a highly radiopaque material for clear fluoroscopic visualization





URETERAL STENTS



| Product Description | Part Number | |
|--------------------------|-------------|----------------|
| 4.9Fr | | Packaged 1/box |
| 4.9Fr x 14cm | CAS4914 | |
| 4.9Fr x 20cm | CAS4920 | |
| 4.9Fr x 22cm | CAS4922 | |
| 4.9Fr x 24cm | CAS4924 | |
| 4.9Fr x 26cm | CAS4926 | |
| 4.9Fr x 28cm | CAS4928 | |
| 4.9Fr x 22-32cm Variable | CAS49VR | |
| <hr/> | | |
| 6Fr | | Packaged 1/box |
| 6Fr x 20cm | CAS6020 | |
| 6Fr x 22cm | CAS6022 | |
| 6Fr x 24cm | CAS6024 | |
| 6Fr x 26cm | CAS6026 | |
| 6Fr x 28cm | CAS6028 | |
| 6Fr x 22-32cm Variable | CAS60VR | |
| <hr/> | | |
| 7Fr | | Packaged 1/box |
| 7Fr x 20cm | CAS7020 | |
| 7Fr x 22cm | CAS7022 | |
| 7Fr x 24cm | CAS7024 | |
| 7Fr x 26cm | CAS7026 | |
| 7Fr x 28cm | CAS7028 | |
| 7Fr x 22-32cm Variable | CAS70VR | |
| <hr/> | | |
| 8Fr | | Packaged 1/box |
| 8Fr x 20cm | CAS8020 | |
| 8Fr x 22cm | CAS8022 | |
| 8Fr x 24cm | CAS8024 | |
| 8Fr x 26cm | CAS8026 | |
| 8Fr x 28cm | CAS8028 | |
| 8Fr x 22-32cm Variable | CAS80VR | |

www.dornier.com

1-800-DORNIER (367-6437)
dornierdirect@dornier.com

© 2020 DORNIER MEDTECH.
ALL RIGHTS RESERVED. SUBJECT TO
CHANGE WITHOUT NOTICE.

DMTA668-042020 REV A EN